

Cleaners and Disinfectants for SonoSite Systems and Transducers

Table of Contents

Introduction	1
Systems	2
Standard Transducers, A-H	10
Standard Transducers, I-Z	17
Specialty Transducers	24
Cleaner and Disinfectant Information	31

Introduction

SonoSite® tests cleaners and disinfectants for use with its systems and transducers. Results are in the tables below. See the ultrasound system user guide for cleaning and disinfecting instructions and for additional warnings and cautions.

See “[Cleaner and Disinfectant Information](#)” for country of origin, type, and active ingredient.

For cleaners and disinfectants compatible with accessories, see Cleaners and Disinfectants for SonoSite Accessories.

For cleaners and disinfectants compatible with TEE transducers, see the TEE Care Card.

For cleaners and disinfectants compatible with surgical transducers, see the Surgical Transducer Guide.

WARNING: Cleaners and disinfectants listed are recommended by SonoSite for safe use with SonoSite products. Refer to the disinfectant label instructions for guidance on disinfection effectiveness and appropriate clinical use.

This document does not have the following regulatory information for disinfectants:

- EPA registration
- FDA 510(k) clearance (liquid sterilant, high level disinfectant)
- CE approval

Before using a disinfectant, confirm that its regulatory status is appropriate for your jurisdiction and use. Verify expiration dates on chemicals.

When disposing of chemicals, follow manufacturer recommendations and EPA regulations.

Systems

The exterior surface of the ultrasound system can be cleaned and disinfected using a recommended cleaner or disinfectant.

Table 1: Cleaner and Disinfectant Compatibility with Sonosite Systems

Cleaner or Disinfectant	✓ - Acceptable for use with system ⊗ - Not Acceptable (Do not use) — Untested (Do not use)						
	NanoMaxx™	S Series™	M-Turbo® M-OB/GYN Office™	MicroMaxx®	TITAN®	SonoSite® 180PLUS™	iLook® 15 iLook® 25
AbcoCide 14	—	—	—	✓	✓	⊗	✓
Accel Plus	—	—	—	⊗	⊗	⊗	⊗
Accel TB Wipes	—	—	—	⊗	⊗	⊗	⊗
Accel Wipes	—	—	—	✓	✓	⊗	⊗
Aidal Plus	—	—	—	⊗	⊗	⊗	⊗
Airkem A-33®	✓	✓	✓	✓	⊗	⊗	✓
Alkacide	—	—	—	✓	✓	✓	⊗
Alkazyme	—	—	—	✓	✓	✓	⊗
Anioxy-Twin	—	—	—	⊗	⊗	⊗	⊗
Anioxyde 1000	⊗	⊗	⊗	✓	⊗	⊗	✓
Aquatabs (1000 PPM)	—	—	—	✓	✓	✓	⊗
Aquatabs (2000 PPM)	—	—	—	⊗	⊗	⊗	⊗
Ascend	—	—	—	✓	✓	✓	✓
Asepti-HB	✓	—	—	⊗	⊗	⊗	⊗

Table 1: Cleaner and Disinfectant Compatibility with Sonosite Systems (continued)

✓ - Acceptable for use with system ⊘ - Not Acceptable (Do not use) — Untested (Do not use)							
Cleaner or Disinfectant	NanoMaxx™	S Series™	M-Turbo® M-OB/GYN Office™	MicroMaxx®	TITAN®	SonoSite® 180PLUS™	iLook® 15 iLook® 25
Asepti-Steryl®	—	⊘	✓	⊘	✓	⊘	⊘
Asepti-Wipes	—	✓	✓	✓	✓	⊘	⊘
Bacillocid rasant	—	—	—	⊘	⊘	✓	✓
Banicide	—	—	—	⊘	⊘	⊘	✓
Bacoban	—	⊘	⊘	⊘	⊘	⊘	⊘
Bacoban WB	—	✓	✓	⊘	⊘	✓	⊘
Betadine	—	—	—	—	—	—	—
Bleach	—	—	—	⊘	⊘	⊘	⊘
Cavicide	—	—	—	✓	✓	⊘	⊘
CaviWipes™	⊘	⊘	⊘	✓	✓	✓	⊘
Chlor-Clean	—	—	—	⊘	⊘	⊘	⊘
Cidalkan®	⊘	⊘	✓	✓	✓	⊘	⊘
Cidalkan® Lingettes	—	—	—	✓	✓	✓	⊘
Cidex®	—	✓	✓	✓	✓	✓	✓
Cidex® OPA	✓	—	—	✓	✓	⊘	✓
Cidex® Plus	—	✓	✓	✓	✓	✓	✓

Table 1: Cleaner and Disinfectant Compatibility with Sonosite Systems (continued)

✓ - Acceptable for use with system ⊘ - Not Acceptable (Do not use) — Untested (Do not use)							
Cleaner or Disinfectant	NanoMaxx™	S Series™	M-Turbo® M-OB/GYN Office™	MicroMaxx®	TITAN®	SonoSite® 180PLUS™	iLook® 15 iLook® 25
Cleanisept Wipe	—	✓	✓	⊘	⊘	⊘	✓
Clorox Wipes	—	—	—	⊘	⊘	⊘	⊘
Control III	—	—	—	✓	✓	⊘	⊘
Coverage™ Spray	—	⊘	⊘	⊘	✓	⊘	⊘
DentaSept	—	—	—	⊘	⊘	✓	⊘
DisCide® Ultra Disinfecting Tow- elettes	—	—	⊘	⊘	⊘	⊘	⊘
DisOPA	—	—	✓	⊘	⊘	⊘	✓
Dispatch	—	—	—	⊘	⊘	⊘	⊘
Dynacide PA	—	—	—	✓	✓	✓	✓
Echo Clean Lingettes	⊘	—	—	—	—	—	—
End-Bac® II	—	⊘	✓	✓	✓	⊘	⊘
Endosporine	—	⊘	✓	✓	⊘	⊘	✓
EndoZime AW® Plus	—	—	—	✓	✓	✓	✓
Envirocide®	—	—	—	✓	✓	⊘	⊘
Enzol	—	—	—	⊘	✓	✓	✓

Table 1: Cleaner and Disinfectant Compatibility with Sonosite Systems (continued)

✓ - Acceptable for use with system ⊘ - Not Acceptable (Do not use) — Untested (Do not use)							
Cleaner or Disinfectant	NanoMaxx™	S Series™	M-Turbo® M-OB/GYN Office™	MicroMaxx®	TITAN®	SonoSite® 180PLUS™	iLook® 15 iLook® 25
Ethyl Alcohol	⊘	—	—	⊘	⊘	⊘	⊘
Expose™ II	—	—	—	✓	✓	✓	⊘
Gigasept AF	—	—	—	⊘	⊘	⊘	⊘
Gigasept FF	—	—	—	⊘	⊘	⊘	⊘
Glutaraldehyde SDS	—	—	—	✓	✓	⊘	⊘
Hexanios	—	—	—	✓	✓	✓	✓
Hi Tor Plus	—	—	—	⊘	⊘	✓	⊘
Hibiclens	—	—	—	✓	✓	✓	✓
Hydrogen Peroxide (3%)	—	—	✓	⊘	⊘	⊘	⊘
Isopropyl Alcohol (100%)	⊘	—	—	⊘	⊘	⊘	⊘
Isopropyl Alcohol (70%)	—	—	—	—	—	—	—
Kodan® Tücher	—	—	⊘	⊘	⊘	⊘	⊘
Kohrsolin® FF	—	—	—	⊘	⊘	⊘	⊘
Korsolex® Basic	—	—	—	✓	✓	⊘	⊘
Korsolex® Extra	—	—	—	⊘	⊘	✓	✓
Korsolex® FF	—	—	—	✓	—	—	—

Table 1: Cleaner and Disinfectant Compatibility with Sonosite Systems (continued)

✓ - Acceptable for use with system ⊘ - Not Acceptable (Do not use) — Untested (Do not use)							
Cleaner or Disinfectant	NanoMaxx™	S Series™	M-Turbo® M-OB/GYN Office™	MicroMaxx®	TITAN®	SonoSite® 180PLUS™	iLook® 15 iLook® 25
Lem-O-Quat	—	—	—	⊘	⊘	⊘	⊘
LpHse	—	—	—	✓	✓	⊘	⊘
Lysol®	—	⊘	—	⊘	⊘	⊘	⊘
Lysol® IC Spray	—	—	—	✓	✓	⊘	⊘
MadaCide 1	⊘	✓	⊘	⊘	⊘	⊘	⊘
Matar	—	—	—	⊘	⊘	⊘	⊘
MetriCide 14	—	—	—	✓	✓	✓	✓
MetriCide® 28	—	—	—	✓	✓	✓	✓
MetriCide OPA Plus	✓	⊘	✓	⊘	⊘	⊘	✓
MetriZyme	—	—	—	✓	✓	✓	✓
Mikrobac forte	—	—	—	✓	✓	✓	✓
Mikrozid Wipes	—	—	⊘	✓	✓	⊘	⊘
Nuclean	—	—	—	⊘	⊘	✓	✓
PerCept RTU Wipes™	✓	—	✓	✓	✓	⊘	⊘
Precise	—	—	—	⊘	⊘	⊘	⊘
Prevention	—	—	—	⊘	⊘	⊘	⊘

Table 1: Cleaner and Disinfectant Compatibility with Sonosite Systems (continued)

✓ - Acceptable for use with system ⊘ - Not Acceptable (Do not use) — Untested (Do not use)							
Cleaner or Disinfectant	NanoMaxx™	S Series™	M-Turbo® M-OB/GYN Office™	MicroMaxx®	TITAN®	SonoSite® 180PLUS™	iLook® 15 iLook® 25
Rely+On™ PeraSafe™	✓	⊘	⊘	⊘	✓	⊘	✓
Ruthless®	—	—	—	✓	✓	⊘	⊘
Sagrosept Wipes	—	—	—	⊘	⊘	⊘	⊘
Salvanios pH 7	—	—	—	✓	✓	✓	✓
Sani-Cloth® HB	—	✓	✓	⊘	⊘	⊘	⊘
Sani-Cloth® Plus	⊘	✓	✓	⊘	⊘	⊘	✓
Sekusept® aktiv	—	—	—	⊘	✓	⊘	⊘
Sklar	—	—	—	⊘	⊘	⊘	⊘
Sporicidin®	⊘	⊘	⊘	⊘	⊘	✓	⊘
Sporicidin® Wipes	—	—	—	⊘	⊘	⊘	⊘
Staphene®	—	—	⊘	⊘	✓	⊘	⊘
Steranios	—	—	—	✓	✓	⊘	✓
Steranios 20%	✓	⊘	⊘	⊘	✓	⊘	✓
Super Sani-Cloth®	—	⊘	⊘	⊘	⊘	⊘	⊘
T-Spray™	⊘	⊘	⊘	⊘	⊘	⊘	⊘
T-Spray II™	✓	✓	✓	✓	✓	✓	✓

Table 1: Cleaner and Disinfectant Compatibility with Sonosite Systems (continued)

✓ - Acceptable for use with system ⊘ - Not Acceptable (Do not use) — Untested (Do not use)							
Cleaner or Disinfectant	NanoMaxx™	S Series™	M-Turbo® M-OB/GYN Office™	MicroMaxx®	TITAN®	SonoSite® 180PLUS™	iLook® 15 iLook® 25
TASK 105	—	—	—	✓	✓	✓	✓
TBQ	—	—	—	⊘	⊘	⊘	⊘
Theracide Plus	—	—	—	⊘	✓	⊘	⊘
Theracide Plus Wipes	—	✓	✓	✓	✓	⊘	⊘
Tor	—	—	—	⊘	⊘	⊘	⊘
Transeptic®	—	—	—	⊘	✓	⊘	⊘
Tristel	—	—	—	✓	✓	⊘	⊘
Tristel Duo	—	—	—	—	—	—	—
Tristel Solo	✓	✓	⊘	✓	⊘	✓	✓
Tristel Wipes	—	⊘	✓	✓	✓	⊘	⊘
Vesphene IIse	—	—	—	✓	✓	⊘	⊘
Virex® II 256	✓	—	—	⊘	⊘	✓	⊘
Virex® TB	⊘	⊘	⊘	⊘	⊘	⊘	⊘
Virox® 5 RTU Wipes	⊘	✓	✓	✓	✓	✓	⊘
Virufen	—	—	—	✓	✓	✓	✓
Wavicide -01	—	—	—	⊘	✓	⊘	✓

Table 1: Cleaner and Disinfectant Compatibility with Sonosite Systems (continued)

✓ - Acceptable for use with system ⊘ - Not Acceptable (Do not use) — Untested (Do not use)							
Cleaner or Disinfectant	NanoMaxx™	S Series™	M-Turbo® M-OB/GYN Office™	MicroMaxx®	TITAN®	SonoSite® 180PLUS™	iLook® 15 iLook® 25
Wavicide -06	—	—	—	⊘	⊘	✓	⊘
Wet Wipe Disinfection	—	—	—	✓	✓	✓	✓
Wex-Cide	—	—	—	✓	✓	✓	✓

Standard Transducers, A-H

The transducer, scanning surface, and cable can be cleaned and disinfected using a recommended cleaner or disinfectant.

Table 2: Cleaner and Disinfectant Compatibility with SonoSite Standard Transducers (A-H)

✓ - Acceptable for use with system ⊗ - Not Acceptable (Do not use) — Untested (Do not use)									
Cleaner or Disinfectant	C11 C11e C11x	C11n	C11 Cast Urethane ¹	C15 C15e	C15 Cast Urethane ²	C60 C60e C60x	C60n	D2 D2x	HFL38 HFL38x
AbcoCide 14	✓	—	✓	✓	✓	✓	—	—	✓
Accel Plus	⊗	—	⊗	⊗	⊗	⊗	—	—	⊗
Accel TB Wipes	⊗	—	⊗	⊗	⊗	⊗	—	—	⊗
Accel Wipes	✓	—	✓	✓	✓	✓	—	—	✓
Aidal Plus	✓	—	✓	✓	✓	✓	—	—	✓
Airkem A-33®	✓	✓	✓	✓	✓	✓	✓	✓	✓
Alkacide	✓	—	✓	✓	✓	✓	—	—	✓
Alkazyme	✓	—	✓	✓	✓	✓	—	—	✓
Anioxy-Twin	⊗	—	⊗	⊗	⊗	⊗	—	—	⊗
Anioxyde 1000	⊗	⊗	⊗	⊗	⊗	⊗	⊗	✓	⊗
Aquatabs (1000)	✓	—	✓	✓	✓	✓	—	—	⊗
Aquatabs (2000)	✓	—	✓	✓	✓	✓	—	—	⊗
Aquatabs (5000)	⊗	—	⊗	⊗	⊗	⊗	—	—	⊗
Ascend	✓	—	✓	✓	✓	✓	—	—	✓
Asepti-HB	✓	✓	✓	✓	✓	✓	✓	—	✓
Asepti-Steryl®	✓	—	⊗	✓	⊗	✓	—	—	✓
Asepti-Wipes	✓	—	✓	✓	✓	✓	—	—	✓

Table 2: Cleaner and Disinfectant Compatibility with SonoSite Standard Transducers (A-H) (continued)

✓ - Acceptable for use with system ⊗ - Not Acceptable (Do not use) — Untested (Do not use)									
Cleaner or Disinfectant	C11 C11e C11x	C11n	C11 Cast Urethane ¹	C15 C15e	C15 Cast Urethane ²	C60 C60e C60x	C60n	D2 D2x	HFL38 HFL38x
Bacillocid rasant	✓	—	✓	✓	✓	✓	—	—	✓
Banicide	✓	—	✓	✓	✓	✓	—	—	—
Bacoban	✓	—	✓	✓	✓	✓	—	✓	✓
Bacoban WB	✓	—	✓	✓	✓	✓	—	✓	✓
Betadine	—	—	—	—	—	—	—	—	—
Bleach	✓	—	✓	✓	✓	✓	—	—	✓
Cavicide	✓	—	⊗	✓	⊗	✓	—	—	✓
CaviWipes™	✓	✓	✓	✓	✓	✓	✓	⊗	✓
Chlor-Clean	✓	—	✓	✓	✓	✓	—	—	⊗
Cidalkan®	✓	—	✓	✓	✓	✓	—	⊗	✓
Cidalkan® Lingettes	✓	—	✓	✓	✓	✓	—	⊗	✓
Cidex®	✓	—	✓	✓	✓	✓	—	—	✓
Cidex® OPA	✓	✓ ³	✓	✓	✓	✓	✓ ³	✓	✓
Cidex® Plus	✓	—	✓	✓	✓	✓	—	—	✓
Cleanisept Wipe	✓	—	✓	✓	✓	✓	—	✓	✓
Clorox Wipes	✓	—	✓	✓	✓	✓	—	—	✓
Control III	⊗	—	✓	✓	✓	✓	—	—	✓
Coverage™ Spray	⊗	⊗	✓	✓	✓	✓	✓	—	✓

Table 2: Cleaner and Disinfectant Compatibility with SonoSite Standard Transducers (A-H) (continued)

✓ - Acceptable for use with system ⊗ - Not Acceptable (Do not use) — Untested (Do not use)									
Cleaner or Disinfectant	C11 C11e C11x	C11n	C11 Cast Urethane ¹	C15 C15e	C15 Cast Urethane ²	C60 C60e C60x	C60n	D2 D2x	HFL38 HFL38x
DentaSept	⊗	—	⊗	⊗	⊗	⊗	—	—	⊗
DisCide® Ultra Disinfecting Tow- elettes	✓	—	⊗	✓	⊗	✓	—	—	✓
DisOPA	✓	—	✓	✓	✓	✓	—	✓	✓
Dispatch	✓	—	✓	✓	✓	✓	—	✓	✓
Dynacide PA	✓	—	✓	✓	✓	✓	—	✓	✓
Echo Clean Lingettes	—	⊗	—	—	—	—	⊗	—	—
End-Bac® II	✓	—	⊗	✓	⊗	✓	—	—	✓
Endosporine	✓	✓	✓	✓	✓	✓	✓	✓	✓
EndoZime AW® Plus	✓	—	✓	✓	✓	✓	—	—	✓
Envirocide®	⊗	⊗	✓	✓	✓	✓	✓	—	—
Enzol	✓	—	✓	✓	✓	✓	—	—	✓
Ethyl Alcohol	⊗	✓	⊗	⊗	⊗	⊗	✓	—	⊗
Expose™	✓	—	✓	✓	✓	✓	—	—	✓
Gigasept AF	✓	—	✓	✓	✓	✓	—	—	✓
Gigasept FF	⊗	—	⊗	⊗	⊗	⊗	—	—	⊗
Glutaraldehyde SDS	✓	—	✓	✓	✓	✓	—	—	—
Hexanios	✓	—	✓	✓	✓	✓	—	—	✓
Hi Tor Plus	⊗	—	✓	✓	✓	✓	—	—	✓

Table 2: Cleaner and Disinfectant Compatibility with SonoSite Standard Transducers (A-H) (continued)

✓ - Acceptable for use with system ⊗ - Not Acceptable (Do not use) — Untested (Do not use)									
Cleaner or Disinfectant	C11 C11e C11x	C11n	C11 Cast Urethane ¹	C15 C15e	C15 Cast Urethane ²	C60 C60e C60x	C60n	D2 D2x	HFL38 HFL38x
Hibiclens	✓	—	✓	✓	✓	✓	—	—	✓
Hydrogen Peroxide (3%)	✓	✓	✓	✓	✓	✓	✓	✓	✓
Isopropyl Alcohol (100%)	⊗	⊗	⊗	⊗	⊗	⊗	⊗	—	⊗
Isopropyl Alcohol (70%)	—	—	—	—	—	—	—	—	—
Kodan® Tücher	✓	—	⊗	✓	⊗	✓	—	—	✓
Kohrsolin® FF	✓	—	✓	✓	✓	✓	—	—	—
Korsolex® Basic	⊗	—	⊗	⊗	⊗	⊗	—	—	⊗
Korsolex® Extra	✓	—	✓	✓	✓	✓	—	—	✓
Korsolex® FF	—	—	—	—	—	—	—	—	—
Lem-O-Quat	⊗	—	⊗	⊗	⊗	✓	—	✓	⊗
LpHse	✓	—	✓	✓	✓	✓	—	—	✓
Lysol®	⊗	—	⊗	⊗	⊗	⊗	—	—	⊗
Lysol® IC Spray	✓	—	⊗	✓	⊗	✓	—	—	⊗
Madacide 1	⊗	⊗	⊗	✓	⊗	✓	✓	⊗	✓
Matar	✓	—	✓	✓	✓	✓	—	—	—
MetriCide 14	✓	—	✓	✓	✓	✓	—	—	✓
MetriCide® 28	✓	—	✓	✓	✓	✓	—	—	✓
Metricide OPA Plus	✓	✓	✓	✓	✓	✓	✓	✓	✓

Table 2: Cleaner and Disinfectant Compatibility with SonoSite Standard Transducers (A-H) (continued)

✓ - Acceptable for use with system ⊗ - Not Acceptable (Do not use) — Untested (Do not use)									
Cleaner or Disinfectant	C11 C11e C11x	C11n	C11 Cast Urethane ¹	C15 C15e	C15 Cast Urethane ²	C60 C60e C60x	C60n	D2 D2x	HFL38 HFL38x
MetriZyme	✓	—	✓	✓	✓	✓	—	—	✓
Mikrobac forte	✓	—	✓	✓	✓	✓	—	—	✓
Mikrozid Wipes	✓	—	⊗	✓	⊗	✓	—	—	✓
Nuclear	✓	—	✓	✓	✓	✓	—	—	✓
PerCept RTU Wipes™	✓	✓	✓	✓	✓	✓	✓	✓	✓
Precise	⊗	—	⊗	⊗	⊗	⊗	—	—	⊗
Prevention	⊗	—	⊗	⊗	⊗	⊗	—	—	⊗
Rely+On™ PeraSafe™	⊗	⊗	⊗	⊗	⊗	⊗	⊗	✓	⊗
Ruthless®	⊗	—	⊗	✓	⊗	✓	—	—	✓
Sagrosept Wipe	✓	—	⊗	✓	⊗	✓	—	—	✓
Salvanios pH 7	✓	—	✓	✓	✓	✓	—	—	✓
Sani-Cloth® HB	⊗	⊗	✓	✓	✓	✓	✓	—	✓
Sani-Cloth® Plus	✓	⊗	✓	✓	✓	✓	⊗	—	✓
Sekusept®	✓	—	✓	✓	✓	✓	—	✓	✓
Sklar	⊗	—	✓	✓	✓	✓	—	—	✓
Sporicidin®	✓	⊗	✓	✓	✓	✓	⊗	⊗	✓
Sporicidin® Wipes	✓	—	✓	✓	✓	✓	—	—	✓
Staphene®	✓	—	⊗	✓	⊗	✓	—	—	⊗

Table 2: Cleaner and Disinfectant Compatibility with SonoSite Standard Transducers (A-H) (continued)

✓ - Acceptable for use with system ⊗ - Not Acceptable (Do not use) — Untested (Do not use)									
Cleaner or Disinfectant	C11 C11e C11x	C11n	C11 Cast Urethane ¹	C15 C15e	C15 Cast Urethane ²	C60 C60e C60x	C60n	D2 D2x	HFL38 HFL38x
Steranios	✓		✓	✓	✓	✓		—	✓
Steranios 20%	⊗	⊗	⊗	⊗	⊗	⊗	⊗	✓	⊗
Super Sani-Cloth®	⊗	⊗	⊗	⊗	⊗	⊗	⊗	—	⊗
T-Spray™	⊗	⊗	✓	✓	✓	✓	✓	—	✓
T-Spray II™	✓	✓ ³	✓	✓	✓	✓	✓ ³	✓	✓
TASK 105	✓	—	✓	✓	✓	✓	—	—	✓
TBQ	✓	—	✓	✓	✓	✓	—	—	✓
Theracide Plus	✓	—	✓	✓	✓	✓	—	—	✓
Theracide Plus Wipes	✓	—	✓	✓	✓	✓	—	—	✓
Tor	⊗	—	✓	✓	✓	✓	—	—	✓
Transeptic®	⊗	—	⊗	⊗	⊗	⊗	—	—	⊗
Tristel	✓	—	⊗	✓	⊗	✓	—	✓	✓
Tristel Duo	⊗	—	⊗	⊗	⊗	⊗	—	—	⊗
Tristel Solo	✓	✓	✓	✓	✓	✓	✓	✓	✓
Tristel Wipes	⊗	—	⊗	⊗	⊗	⊗	—	✓	⊗
Vesphene IIse	✓	—	✓	✓	✓	✓	—	—	✓
Virex® II 256	✓	✓	✓	✓	✓	✓	✓	—	✓
Virex® TB	⊗	⊗	✓	✓	✓	✓	✓	—	✓

Table 2: Cleaner and Disinfectant Compatibility with SonoSite Standard Transducers (A-H) (continued)

✓ - Acceptable for use with system ⊗ - Not Acceptable (Do not use) — Untested (Do not use)									
Cleaner or Disinfectant	C11 C11e C11x	C11n	C11 Cast Urethane ¹	C15 C15e	C15 Cast Urethane ²	C60 C60e C60x	C60n	D2 D2x	HFL38 HFL38x
Virox® 5 RTU Wipes	✓	✓	✓	✓	✓	✓	✓	✓	✓
Virufen	✓	—	✓	✓	✓	✓	—	✓	✓
Wavicide -01	⊗	—	⊗	⊗	⊗	⊗	—	—	⊗
Wavicide -06	✓	—	⊗	✓	⊗	✓	—	—	✓
Wet Wipe Disinfection	✓	—	✓	✓	✓	✓	—	✓	✓
Wex-Cide	✓	—	✓	✓	✓	✓	—	—	✓

1 Part number: P01960-xx, P02645-02, P02645-03, or P03361-10

2 Part number: P00798-03

3 Cleaner or disinfectant has been tested by SonoSite for disinfection effectiveness on SonoSite transducers.

Standard Transducers, I-Z

The transducer, scanning surface, and cable can be cleaned and disinfected using a recommended cleaner or disinfectant.

Table 3: Cleaner and Disinfectant Compatibility with SonoSite Standard Transducers (I-Z)

Cleaner or Disinfectant	L25 L25e L25x	L25n	L38 L38e L38x	L38n	L52 L52e Vet Only	L52x Vet Only	P10 P10x	P17	P21x	P21n
AbcoCide 14	✓	—	✓	—	✓	—	✓	✓	✓	—
Accel Plus	⊘	—	⊘	—	⊘	—	⊘	⊘	⊘	—
Accel TB Wipes	⊘	—	⊘	—	⊘	—	⊘	⊘	⊘	—
Accel Wipes	✓	—	✓	—	✓	—	✓	✓	✓	—
Aidal Plus	✓	—	✓	—	✓	—	✓	✓	✓	—
Airkem A-33®	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Alkacide	✓	—	✓	—	✓	—	✓	✓	✓	—
Alkazyme	✓	—	✓	—	✓	—	✓	✓	✓	—
Anioxy-Twin	⊘	—	⊘	—	⊘	—	⊘	⊘	⊘	—
Anioxyde 1000	⊘	⊘	⊘	⊘	✓	⊘	⊘	⊘	⊘	⊘
Aquatabs (1000)	✓	—	✓	—	✓	—	✓	✓	✓	—
Aquatabs (2000)	✓	—	✓	—	✓	—	✓	✓	✓	—
Aquatabs (5000)	⊘	—	⊘	—	⊘	—	⊘	⊘	⊘	—
Ascend	✓	—	✓	—	✓	—	✓	✓	✓	—
Asepti-HB	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Asepti-Steryl®	✓	—	✓	—	✓	—	✓	✓	✓	—

Table 3: Cleaner and Disinfectant Compatibility with SonoSite Standard Transducers (I-Z) (continued)

✓ - Acceptable for use with system ⊗ - Not Acceptable (Do not use) — Untested (Do not use)										
Cleaner or Disinfectant	L25 L25e L25x	L25n	L38 L38e L38x	L38n	L52 L52e Vet Only	L52x Vet Only	P10 P10x	P17	P21x	P21n
Asepti-Wipes	✓	—	✓	—	✓	—	✓	✓	✓	—
Bacillocid rasant	✓	—	✓	—	✓	—	✓	✓	✓	—
Banicide	✓	—	✓	—	✓	—	✓	✓	✓	—
Bacoban	✓	—	✓	—	✓	—	✓	✓	✓	—
Bacoban WB	✓	—	✓	—	✓	—	✓	✓	✓	—
Betadine	—	—	—	—	—	—	—	—	—	—
Bleach	✓	—	✓	—	✓	—	✓	✓	✓	—
Cavicide	✓	—	✓	—	✓	—	✓	✓	✓	—
CaviWipes™	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Chlor-Clean	✓	—	✓	—	✓	—	✓	✓	✓	—
Cidalkan®	✓	—	✓	—	✓	—	✓	✓	✓	—
Cidalkan® Lingettes	✓	—	✓	—	✓	—	✓	✓	✓	—
Cidex®	✓	—	✓	—	✓	—	✓	✓	✓	—
Cidex® OPA	✓	✓ ¹	✓	✓ ¹	✓	✓	✓	✓	✓	✓ ¹
Cidex® Plus	✓	—	✓	—	✓	—	✓	✓	✓	—
Cleanisept Wipe	✓	—	✓	—	✓	—	✓	✓	✓	—
Clorox Wipes	✓	—	✓	—	✓	—	✓	✓	✓	—
Control III	⊗	—	✓	—	✓	—	✓	✓	✓	—

Table 3: Cleaner and Disinfectant Compatibility with SonoSite Standard Transducers (I-Z) (continued)

✓ - Acceptable for use with system ⊗ - Not Acceptable (Do not use) — Untested (Do not use)										
Cleaner or Disinfectant	L25 L25e L25x	L25n	L38 L38e L38x	L38n	L52 L52e Vet Only	L52x Vet Only	P10 P10x	P17	P21x	P21n
Coverage™ Spray	⊗	⊗	✓	✓	✓	—	✓	✓	✓	✓
DentaSept	⊗	—	⊗	—	✓	—	⊗	⊗	⊗	—
DisCide® Ultra Disinfecting Tow- elettes	✓	—	✓	—	✓	—	✓	✓	✓	—
DisOPA	✓	—	✓	—	✓	—	✓	✓	✓	—
Dispatch	✓	—	✓	—	✓	—	✓	✓	✓	—
Dynacide PA	✓	—	✓	—	✓	—	✓	✓	✓	—
Echo Clean Lingettes	—	⊗	—	⊗	—	—	—	—	—	⊗
End-Bac® II	✓	—	✓	—	✓	—	✓	✓	✓	—
Endosporine	✓	✓	✓	✓	✓	⊗	✓	✓	✓	✓
EndoZime AW® Plus	✓	—	✓	—	✓	—	✓	✓	✓	—
Envirocide®	⊗	⊗	✓	✓	✓	⊗	✓	✓	✓	✓
Enzol	✓	—	✓	—	✓	—	✓	✓	✓	—
Ethyl Alcohol	⊗	✓	⊗	✓	⊗	✓	⊗	⊗	⊗	✓
Expose™	✓	—	✓	—	✓	—	✓	✓	✓	—
Gigasept AF	✓	—	✓	—	✓	—	✓	✓	✓	—
Gigasept FF	⊗	—	⊗	—	✓	—	⊗	⊗	⊗	—
Glutaraldehyde SDS	✓	—	✓	—	✓	—	✓	✓	✓	—

Table 3: Cleaner and Disinfectant Compatibility with SonoSite Standard Transducers (I-Z) (continued)

✓ - Acceptable for use with system ⊗ - Not Acceptable (Do not use) — Untested (Do not use)										
Cleaner or Disinfectant	L25 L25e L25x	L25n	L38 L38e L38x	L38n	L52 L52e L52x Vet Only	L52x Vet Only	P10 P10x	P17	P21x	P21n
Hexanios	✓	—	✓	—	✓	—	✓	✓	✓	—
Hi Tor Plus	⊗	—	✓	—	✓	—	✓	✓	✓	—
Hibiclens	✓	—	✓	—	✓	—	✓	✓	✓	—
Hydrogen Peroxide (3%)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Isopropyl Alcohol (100%)	⊗	⊗	⊗	⊗	⊗	⊗	⊗	⊗	⊗	⊗
Isopropyl Alcohol (70%)	—	—	—	—	—	—	—	—	—	—
Kodan® Tücher	✓	—	✓	—	✓	—	✓	✓	✓	—
Kohrsolin® FF	✓	—	✓	—	✓	—	✓	✓	✓	—
Korsolex® Basic	⊗	—	⊗	—	✓	—	⊗	⊗	⊗	—
Korsolex® Extra	✓	—	✓	—	✓	—	✓	✓	✓	—
Korsolex® FF	—	—	—	—	—	—	—	—	—	—
Lem-O-Quat	⊗	—	⊗	—	⊗	—	⊗	⊗	⊗	—
LpHse	✓	—	✓	—	✓	—	✓	✓	✓	—
Lysol®	⊗	—	⊗	—	⊗	—	⊗	⊗	⊗	—
Lysol® IC Spray	✓	—	✓	—	✓	—	✓	✓	✓	—
Madacide 1	⊗	⊗	✓	✓	✓	✓	✓	✓	✓	✓
Matar	✓	—	✓	—	✓	—	✓	✓	✓	—
MetriCide 14	✓	—	✓	—	✓	—	✓	✓	✓	—

Table 3: Cleaner and Disinfectant Compatibility with SonoSite Standard Transducers (I-Z) (continued)

✓ - Acceptable for use with system ⊗ - Not Acceptable (Do not use) — Untested (Do not use)										
Cleaner or Disinfectant	L25 L25e L25x	L25n	L38 L38e L38x	L38n	L52 L52e L52x Vet Only	L52x Vet Only	P10 P10x	P17	P21x	P21n
MetriCide® 28	✓	—	✓	—	✓	—	✓	✓	✓	—
Metricide OPA Plus	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
MetriZyme	✓	—	P	—	✓	—	✓	✓	✓	—
Mikrobac forte	✓	—	✓	—	✓	—	✓	✓	✓	—
Mikrozid Wipes	✓	—	✓	—	✓	—	✓	✓	✓	—
Nuclear	✓	—	✓	—	✓	—	✓	✓	✓	—
PerCept RTU Wipes™	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Precise	⊗	—	⊗	—	⊗	—	⊗	⊗	⊗	—
Prevention	⊗	—	⊗	—	⊗	—	⊗	⊗	⊗	—
Rely+On™ PeraSafe™	⊗	⊗	⊗	⊗	✓	✓	⊗	⊗	⊗	⊗
Ruthless®	⊗	—	✓	—	✓	—	✓	✓	✓	—
Sagrosept Wipe	✓	—	✓	—	✓	—	✓	✓	✓	—
Salvanios pH 7	✓	—	✓	—	✓	—	✓	✓	✓	—
Sani-Cloth® HB	⊗	⊗	✓	✓	✓	✓	✓	✓	✓	✓
Sani-Cloth® Plus	✓	⊗	✓	⊗	✓	✓	✓	✓	✓	⊗
Sekusept®	✓	—	✓	—	✓	—	✓	✓	✓	—
Sklar	⊗	—	✓	—	✓	—	✓	✓	✓	—
Sporicidin®	✓	⊗	✓	⊗	✓	—	✓	✓	✓	⊗

Table 3: Cleaner and Disinfectant Compatibility with SonoSite Standard Transducers (I-Z) (continued)

✓ - Acceptable for use with system ⊘ - Not Acceptable (Do not use) — Untested (Do not use)										
Cleaner or Disinfectant	L25 L25e L25x	L25n	L38 L38e L38x	L38n	L52 L52e Vet Only	L52x Vet Only	P10 P10x	P17	P21x	P21n
Sporicidin® Wipes	✓	—	✓	⊘	✓	—	✓	✓	✓	⊘
Staphene®	✓	—	✓	—	✓	—	✓	✓	✓	—
Steranios	✓	—	✓	—	✓	—	✓	✓	✓	—
Steranios 20%	⊘	⊘	⊘	⊘	✓	⊘	⊘	⊘	⊘	⊘
Super Sani-Cloth®	⊘	⊘	⊘	⊘	✓	⊘	⊘	⊘	⊘	⊘
T-Spray™	⊘	⊘	✓	✓	✓	✓	✓	✓	✓	✓
T-Spray II™	✓	✓ ¹	✓	✓ ¹	✓	✓	✓	✓	✓	✓ ¹
TASK 105	✓	—	✓	—	✓	—	✓	✓	✓	—
TBQ	✓	—	✓	—	✓	—	✓	✓	✓	—
Theracide Plus	✓	—	✓	—	✓	—	✓	✓	✓	—
Theracide Plus Wipes	✓	—	✓	—	⊘	—	✓	✓	✓	—
Tor	⊘	—	✓	—	✓	—	✓	✓	✓	—
Transeptic®	⊘	—	⊘	—	⊘	—	⊘	⊘	⊘	—
Tristel	✓	—	✓	—	✓	—	✓	✓	✓	—
Tristel Duo	⊘	—	⊘	—	⊘	—	⊘	⊘	⊘	—
Tristel Solo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Tristel Wipes	⊘	—	⊘	—	✓	—	⊘	⊘	⊘	—
Vesphene IIse	✓	—	✓	—	✓	—	✓	✓	✓	—

Table 3: Cleaner and Disinfectant Compatibility with SonoSite Standard Transducers (I-Z) (continued)

✓ - Acceptable for use with system ⊗ - Not Acceptable (Do not use) — Untested (Do not use)										
Cleaner or Disinfectant	L25 L25e L25x	L25n	L38 L38e L38x	L38n	L52 L52e Vet Only	L52x Vet Only	P10 P10x	P17	P21x	P21n
Virex® II 256	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Virex® TB	⊗	⊗	✓	✓	✓	✓	✓	✓	✓	✓
Virox® 5 RTU Wipes	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Virufen	✓	—	✓	—	✓	—	✓	✓	✓	—
Wavicide -01	⊗	—	⊗	—	⊗	—	⊗	⊗	⊗	—
Wavicide -06	✓	—	✓	—	✓	—	✓	✓	✓	—
Wet Wipe Disinfection	✓	—	✓	—	✓	—	✓	✓	✓	—
Wex-Cide	✓	—	✓	—	✓	—	✓	✓	✓	—

1 Cleaner or disinfectant has been tested by SonoSite for disinfection effectiveness on SonoSite transducers.

Specialty Transducers

The transducer, scanning surface, and cable can be cleaned and disinfected using a recommended cleaner or disinfectant.

Table 4: Cleaner and Disinfectant Compatibility with SonoSite Specialty Transducers

✓ - Acceptable for use with system ⊗ - Not Acceptable (Do not use) — Untested (Do not use)					
Cleaner or Disinfectant	C8 C8e	HST	ICT ICTx	SLA SLAx	SLT
AbcoCide 14	✓	—	✓	✓	✓
Accel Plus	⊗	—	⊗	⊗	⊗
Accel TB Wipes	⊗	—	⊗	⊗	⊗
Accel Wipes	✓	—	✓	✓	✓
Aidal Plus	✓	✓	✓	✓	✓
Airkem A-33®	✓	—	✓	✓	—
Alkacide	✓	—	✓	✓	✓
Alkazyme	✓	⊗	✓	✓	✓
Anioxy-Twin	⊗	—	⊗	⊗	⊗
Anioxyde 1000	⊗	⊗	⊗	⊗	⊗
Aquatabs (1000)	✓	—	✓	✓	✓
Aquatabs (2000)	✓	—	✓	✓	✓
Aquatabs (5000)	⊗	—	⊗	⊗	⊗
Ascend	✓	—	✓	✓	✓
Asepti-HB	✓	—	✓	✓	✓
Asepti-Steryl®	✓	—	✓	✓	✓
Asepti-Wipes	✓	—	✓	✓	✓

Table 4: Cleaner and Disinfectant Compatibility with SonoSite Specialty Transducers (continued)

✓ - Acceptable for use with system ⊘ - Not Acceptable (Do not use) — Untested (Do not use)					
Cleaner or Disinfectant	C8 C8e	HST	ICT ICTx	SLA SLAx	SLT
Bacillocid rasant	✓	—	✓	✓	✓
Banicide	✓	✓	✓	✓	✓
Bacoban	✓	—	✓	✓	✓
Bacoban WB	✓	—	✓	✓	✓
Betadine	—	—	—	—	—
Bleach	✓	✓	✓	✓	✓
Cavicide	✓	—	✓	✓	✓
CaviWipes™	✓	—	✓	✓	✓
Chlor-Clean	✓	—	✓	✓	✓
Cidalkan®	✓	—	✓	✓	✓
Cidalkan® Lingettes	✓	—	✓	✓	✓
Cidex®	✓	—	✓	✓	✓
Cidex® OPA	✓	✓	✓	✓	✓
Cidex® Plus	✓	—	✓	✓	✓
Cleanisept Wipe	✓	—	✓	✓	✓
Clorox Wipes	✓	—	✓	✓	✓
Control III	✓	—	✓	✓	✓
Coverage™ Spray	✓	—	✓	✓	✓
DentaSept	⊘	—	⊘	⊘	✓

Table 4: Cleaner and Disinfectant Compatibility with SonoSite Specialty Transducers (continued)

✓ - Acceptable for use with system ⊘ - Not Acceptable (Do not use) — Untested (Do not use)					
Cleaner or Disinfectant	C8 C8e	HST	ICT ICTx	SLA SLAx	SLT
DisCide® Ultra Disinfecting Tow-elettes	✓	—	✓	✓	✓
DisOPA	✓	—	✓	✓	✓
Dispatch	✓	✓	✓	✓	✓
Dynacide PA	✓	—	✓	✓	✓
Echo Clean Lingettes	—	—	—	—	—
End-Bac® II	✓	—	✓	✓	✓
Endosporine	✓	—	✓	✓	—
EndoZime AW® Plus	✓	—	✓	✓	✓
Envirocide®	✓	—	✓	✓	✓
Enzol	✓	✓	✓	✓	✓
Ethyl Alcohol	⊘	—	⊘	⊘	⊘
Expose™	✓	—	✓	✓	✓
Gigasept AF	✓	—	✓	✓	✓
Gigasept FF	⊘	✓	⊘	⊘	✓
Glutaraldehyde SDS	✓	—	✓	✓	✓
Hexanios	✓	—	✓	✓	✓
Hi Tor Plus	✓	—	✓	✓	✓
Hibiclens	✓	—	✓	✓	✓

Table 4: Cleaner and Disinfectant Compatibility with SonoSite Specialty Transducers (continued)

✓ - Acceptable for use with system ⊗ - Not Acceptable (Do not use) — Untested (Do not use)					
Cleaner or Disinfectant	C8 C8e	HST	ICT ICTx	SLA SLAx	SLT
Hydrogen Peroxide (3%)	✓	—	✓	✓	✓
Isopropyl Alcohol (100%)	⊗	—	⊗	⊗	⊗
Isopropyl Alcohol (70%)	—	—	—	—	—
Kodan® Tücher	✓	—	✓	✓	✓
Kohrsolin® FF	✓	—	✓	✓	✓
Korsolex® Basic	⊗	—	⊗	⊗	✓
Korsolex® Extra	✓	—	✓	✓	✓
Korsolex® FF	—	—	—	—	—
Lem-O-Quat	⊗	—	⊗	⊗	⊗
LpHse	✓	—	✓	✓	✓
Lysol®	⊗	—	⊗	⊗	⊗
Lysol® IC Spray	✓	—	✓	✓	✓
Madacide 1	✓	—	✓	✓	✓
Matar	✓	—	✓	✓	✓
MetriCide 14	✓	—	✓	✓	✓
MetriCide® 28	✓	✓	✓	✓	✓
Metricide OPA Plus	✓	—	✓	✓	—
MetriZyme	✓	✓	✓	✓	✓
Mikrobac forte	✓	—	✓	✓	✓

Table 4: Cleaner and Disinfectant Compatibility with SonoSite Specialty Transducers (continued)

✓ - Acceptable for use with system ⊘ - Not Acceptable (Do not use) — Untested (Do not use)					
Cleaner or Disinfectant	C8 C8e	HST	ICT ICTx	SLA SLAx	SLT
Mikrozid Wipes	✓	—	✓	✓	✓
Nuclear	✓	—	✓	✓	✓
PerCept RTU Wipes™	—	—	✓	✓	—
Precise	⊘	—	⊘	⊘	⊘
Prevention	⊘	—	⊘	⊘	⊘
Rely+On™ PeraSafe™	⊘	—	⊘	⊘	—
Ruthless®	✓	⊘	✓	✓	✓
Sagrosept Wipe	✓	—	✓	✓	✓
Salvanios pH 7	✓	—	✓	✓	✓
Sani-Cloth® HB	✓	—	✓	✓	✓
Sani-Cloth® Plus	✓	—	✓	✓	✓
Sekusept®	✓	—	✓	✓	✓
Sklar	✓	—	✓	✓	✓
Sporicidin®	✓	—	✓	✓	✓
Sporicidin® Wipes	✓	—	✓	✓	✓
Staphene®	✓	—	✓	✓	✓
Steranios	✓	—	✓	✓	✓
Steranios 20%	⊘	—	⊘	⊘	—
Super Sani-Cloth®	⊘	⊘	⊘	⊘	✓

Table 4: Cleaner and Disinfectant Compatibility with SonoSite Specialty Transducers (continued)

✓ - Acceptable for use with system ⊗ - Not Acceptable (Do not use) — Untested (Do not use)					
Cleaner or Disinfectant	C8 C8e	HST	ICT ICTx	SLA SLAx	SLT
T-Spray™	✓	—	✓	✓	✓
T-Spray II™	✓	✓	✓	✓	✓
TASK 105	✓	—	✓	✓	✓
TBQ	✓	—	✓	✓	✓
Theracide Plus	✓	—	✓	✓	✓
Theracide Plus Wipes	✓	⊗	✓	✓	⊗
Tor	✓	—	✓	✓	✓
Transeptic®	⊗	—	⊗	⊗	⊗
Tristel	✓	—	✓	✓	✓
Tristel Duo	✓	—	⊗	⊗	⊗
Tristel Solo	✓	—	✓	✓	✓
Tristel Wipes	⊗	—	⊗	⊗	✓
Vesphene IIse	✓	—	✓	✓	✓
Virex® II 256	✓	—	✓	✓	✓
Virex® TB	✓	—	✓	✓	✓
Virox® 5 RTU Wipes	✓	—	✓	✓	✓
Virufen	✓	—	✓	✓	✓
Wavicide -01	⊗	✓	⊗	⊗	⊗
Wavicide -06	✓	—	✓	✓	✓

Table 4: Cleaner and Disinfectant Compatibility with SonoSite Specialty Transducers (continued)

✓ - Acceptable for use with system
⊗ - Not Acceptable (Do not use)
— Untested (Do not use)

Cleaner or Disinfectant	C8 C8e	HST	ICT ICTx	SLA SLAx	SLT
Wet Wipe Disinfection	✓	—	✓	✓	✓
Wex-Cide	✓	—	✓	✓	✓

Cleaner and Disinfectant Information

For additional chemical composition information, refer to the cleaner or disinfectant Material Safety Data Sheet (MSDS).

Table 5: Cleaner and Disinfectant Information

Cleaner or Disinfectant	Country of Origin	Type	Active Ingredient
AbcoCide 14	USA	Liquid	Glutaraldehyde
Accel Plus	CAN	Liquid	Hydrogen Peroxide
Accel TB Wipes	CAN	Wipe	Hydrogen Peroxide
Accel Wipes	CAN	Wipe	Hydrogen Peroxide
Aidal Plus	AUS	Liquid	Glutaraldehyde
Airkem A-33®	USA	Liquid	Quaternary Ammonia
Alkacide	FRA	Liquid	Glutaraldehyde
Alkazyme	FRA	Liquid	Quaternary Ammonia
Anioxy-Twin	USA	Liquid	Peracetic Acid
Anioxyde 1000	FRA	Liquid	Peracetic Acid
Aquatabs (1000 PPM)	IRL	Liquid	Chlorine
Aquatabs (2000 PPM)	IRL	Liquid	Chlorine
Ascend	USA	Liquid	Quaternary Ammonia
Asepti-HB	USA	Liquid	Quaternary Ammonia
Asepti-Steryl®	USA	Spray	Ethyl Alcohol
Asepti-Wipes	USA	Wipe	Isopropyl Alcohol
Bacillocid rasant	DEU	Liquid	Glutaraldehyde/Quaternary Ammonia
Banicide	USA	Liquid	Glutaraldehyde
Bacoban	DEU	Liquid	Ethyl Alcohol, Isopropyl Alcohol
Bacoban WB	DEU	Liquid	Benzalkoniumchloride Diethylenglycol

Table 5: Cleaner and Disinfectant Information (continued)

Cleaner or Disinfectant	Country of Origin	Type	Active Ingredient
Betadine	USA	Liquid	Providone-Iodine
Bleach	n/a	Liquid	NaCl Hypochlorite
Cavicide	USA	Liquid	Isopropyl Alcohol
CaviWipes™	USA	Wipes	Isopropyl Alcohol
Chlor-Clean	GBR	Liquid	Sodium Dichloroisocyanurate
Cidalkan®	FRA	Liquid	Alkylamine, Isopropyl Alcohol
Cidalkan® Lingettes	FRA	Wipes	Ethyl Alcohol
Cidex®	USA	Liquid	Glutaraldehyde
Cidex® OPA	USA	Liquid	Ortho-phthaldehyde
Cidex® Plus	USA	Liquid	Glutaraldehyde
Cleanisept Wipe	DEU	Wipes	Quaternary Ammonia
Clorox Wipes	USA	Wipes	Isopropyl Alcohol
Control III	USA	Liquid	Quaternary Ammonia
Coverage™ Spray	USA	Spray	Quaternary Ammonia
DentaSept	FRA	Liquid	Quaternary Ammonia
DisCide® Ultra Disinfecting Tow-elettes	USA	Wipes	Isopropyl Alcohol
DisOPA	JPN	Liquid	Ortho-phthaldehyde
Dispatch	USA	Spray	NaCl Hypochlorite
Dynacide PA	FRA	Liquid	Peracetic Acid
Echo Clean Lingettes	FRA	Wipe	Alkylamine, Isopropyl Alcohol
End-Bac® II	USA	Liquid	Quaternary Ammonia

Table 5: Cleaner and Disinfectant Information (continued)

Cleaner or Disinfectant	Country of Origin	Type	Active Ingredient
Endosporine	FRA	Liquid	Glutaraldehyde
EndoZime AW® Plus	FRA	Liquid	Propanol
Envirocide®	USA	Liquid	Isopropyl Alcohol
Enzol	USA	Cleaner	Ethylene Glycol
Ethyl Alcohol	USA	Liquid	Denatured Ethyl Alcohol 3A
Expose™	USA	Liquid	Isopropyl Alcohol
Gigasept AF	DEU	Liquid	Quaternary Ammonia
Gigasept FF	DEU	Liquid	Bersteinsaure
Glutaraldehyde SDS	USA	Liquid	Glutaraldehyde
Hexanios	FRA	Liquid	Polyhexanide, Quaternary Ammonia
Hi Tor Plus	USA	Liquid	Ammonium Chloride, Quaternary Ammonia, Ethyl Alcohol, Acetic Acid
Hibiclens	USA	Cleaner	Chlorhexidine
Hydrogen Peroxide (3%)	n/a	Liquid	Hydrogen Peroxide
Isopropyl Alcohol (100%)	n/a	Liquid	Isopropyl Alcohol
Isopropyl Alcohol (70%)	n/a	Liquid	Isopropyl Alcohol
Kodan® Tücher	DEU	Liquid	1-Propanol, Isopropyl Alcohol
Kohrsolin® FF	DEU	Liquid	Glutaraldehyde
Korsolex® Basic	DEU	Liquid	Glutaraldehyde
Korsolex® Extra	DEU	Liquid	Isopropyl Alcohol, Quaternary Ammonia, Others
Korsolex® FF	DEU	Liquid	Glutaraldehyde
Lem-O-Quat	USA	Liquid	Quaternary Ammonia

Table 5: Cleaner and Disinfectant Information (continued)

Cleaner or Disinfectant	Country of Origin	Type	Active Ingredient
LpHse	USA	Liquid	Phenolics, Phosphoric Acid, Isopropyl Alcohol, Others
Lysol®	USA	Spray	Ethyl Alcohol
Lysol® IC Spray	USA	Spray	Ethyl Alcohol
Madacide 1	USA	Liquid	Quaternary Ammonia
Matar	USA	Liquid	Phenolics, Ethelene Glycol, Isopropyl Alcohol
MetriCide 14	USA	Liquid	Glutaraldehyde
MetriCide® 28	USA	Liquid	Glutaraldehyde
Metricide OPA Plus	USA	Liquid	ortho-Phthaldehyde
MetriZyme	USA	Cleaner	Propylene Glycol
Mikrobac forte	DEU	Liquid	Ammonium Chloride
Mikrozid Wipes	DEU	Wipe	Ethyl Alcohol, n-Propyl Alcohol
Nuclean	FRA	Spray	Alcohol/Biguanide
PerCept RTU Wipes™	CAN	Wipe	Hydrogen Peroxide
Precise	USA	Spray	Quaternary Ammonia
Prevention	CAN	Liquid	Hydrogen Peroxide
Rely+On™ PeraSafe™	GBR	Liquid	Paracetic Acid
Ruthless®	USA	Spray	Quaternary Ammonia
Sagrosept Wipes	DEU	Wipe	Propanol
Salvanios pH 7	FRA	Liquid	Quaternary Ammonia
Sani-Cloth® HB	USA	Wipe	Quaternary Ammonia
Sani-Cloth® Plus	USA	Wipe	Quaternary Ammonia
Sekusept®	DEU	Liquid	Glutaraldehyde

Table 5: Cleaner and Disinfectant Information (continued)

Cleaner or Disinfectant	Country of Origin	Type	Active Ingredient
Sklar	USA	Liquid	Isopropyl Alcohol
Sporicidin®	USA	Liquid	Phenolics, Glutaraldehyde
Sporicidin® Wipes	USA	Wipe	Phenolics
Staphene®	USA	Spray	Ethyl Alcohol
Steranios	FRA	Liquid	Glutaraldehyde
Steranios 20%	FRA	Liquid	Glutaraldehyde
Super Sani-Cloth®	USA	Wipe	Isopropyl Alcohol
T-Spray™	USA	Spray	Quaternary Ammonia
T-Spray II™	USA	Spray	Alkyl/Chloride
TASK 105	USA	Spray	Quaternary Ammonia
TBQ	USA	Liquid	Alkyl
Theracide Plus	USA	Liquid	Quaternary Ammonia
Theracide Plus Wipes	USA	Wipe	Quaternary Ammonia
Tor	USA	Spray	Quaternary Ammonia
Transeptic®	USA	Cleaner	Alcohol
Tristel	GBR	Liquid	Chlorine Dioxide
Tristel Duo	GBR	Foam	Chlorine Dioxide
Tristel Solo	GBR	Foam	Hexamethylenebiguanide
Tristel Wipes	GBR	Wipe	Chlorine Dioxide
Vesphene IIse	USA	Liquid	Phenolics, Potassium Hydroxide, Sodium Hydroxide
Virex® II 256	USA	Liquid	Ammonium Chloride
Virex® TB	USA	Liquid	Quaternary Ammonia

Table 5: Cleaner and Disinfectant Information (continued)

Cleaner or Disinfectant	Country of Origin	Type	Active Ingredient
Virox® 5 RTU Wipes	CAN	Wipe	Hydrogen Peroxide
Virufen	FRA	Liquid	Quaternary Ammonia
Wavicide -01	USA	Liquid	Glutaraldehyde
Wavicide -06	USA	Liquid	Glutaraldehyde
Wet Wipe Disinfection	DEN	Wipe	Guanidinium-chloride
Wex-Cide	USA	Liquid	Phenolics

iLook, M-OB/GYN Office, M-Turbo, MicroMaxx, NanoMaxx, S Series, SonoSite, SonoSite 180PLUS, TITAN, and the SonoSite logo are either trademarks or registered trademarks of SonoSite, Inc.

SonoSite, Inc.
21919 30th Drive SE
Bothell, WA 98021 USA
T: 1-888-482-9449 or
1-425-951-1200
F: 1-425-951-1201

SonoSite Ltd
Alexander House
40A Wilbury Way
Hitchin, Herts
SG4 0AP UK
T: +44-1462-444800
F: +44-1462-444801

P06703-05 03/2010
Copyright 2009 by
SonoSite, Inc.
All rights reserved.

