






















The following disinfectants are approved for use:

- P.I.-Spray II (replaces T-Spray II)
- Sani Cloth Plus
- Sani Cloth HB
- Protex Wipes
- Cavi Wipes

For Terason systems/transducers:

<p><b>t3000</b></p> 	<p><b>t3200</b></p> 	<p><b>uSmart® 3200T</b></p> 	<p><b>uSmart® 3300</b></p> 
<p><b>12L5</b></p> 	<p><b>15L4</b></p> 	<p><b>15L4</b></p> 	<p><b>15L4</b></p> 
<p><b>5C2</b></p> 	<p><b>6C1</b></p> 	<p><b>5C2</b></p> 	<p><b>5C2</b></p> 
<p><b>4V2</b></p> 	<p><b>8L2</b></p> 	<p><b>4V2</b></p> 	<p><b>4V2</b></p> 
<p><b>12HL7</b></p> 	<p><b>16HL7</b></p> 	<p><b>16HL7</b></p> 	
<p><b>8EC4</b></p> 	<p><b>5C2</b></p> 	<p><b>8EC4</b></p> 